


SXR Series

M - SM - AU Series

B Series



Pediatric | Bicuspids | Bicuspids with deep cervical preps | Bicuspids and small molars | Molars | Molars with deep cervical preps

Slick Bands™ XR Slick Bands™ XR Matrices have non-stick coating and Grap-Tabs™ for easier placement.

thickness 0.0016"						
Height	NA	SXR100 4.6 mm	SXR150 3.8 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm
Refill Quantity		100/pk	60/pk	100/pk	100/pk	60/pk

Slick Bands™ Slick Bands™ Matrices have non-stick coating.

thickness 0.0016"						
Height	NA	SM100 4.6 mm	SM150 3.8 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm
Refill Quantity		100/pk	60/pk	100/pk	100/pk	60/pk

Composi-Tight® Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.

thickness 0.0015							
Height		M050 3.2 mm	M100 4.6 mm	M150 3.8 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm
Refill Quantity		50/pk	100/pk	50/pk	100/pk	100/pk	50/pk

Composi-Tight Gold® The Gold bands are the same as the Composi-Tight® bands.

thickness 0.0015					NA		
Height		AU050 3.2 mm	AU100 4.6 mm	AU150 3.8 mm	NA	AU200 6.4 mm	AU300 6.4 mm
Refill Quantity		50/pk	100/pk	50/pk		100/pk	50/pk

Composi-Tight® Original The original bands have a slightly shorter length.

thickness 0.0015					NA		
Height		B050 3.2 mm	B100 4.6 mm	B150 3.8 mm	NA	B200 6.4 mm	B300 6.4 mm
Refill Quantity		50/pk	100/pk	50/pk		100/pk	50/pk

Composi-Tight® 3D Clear Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.

thickness 0.0020			NA			NA
Height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA
Refill Quantity		100/pk		100/pk	100/pk	