

FOR IMMEDIATE RELEASE

Contact: Jeff Maglin, Director of Marketing

Email: jmaglin@garrisondental.com

Phone: 888-437-0032 x218

Garrison® Dental Solutions Introduces the Next Evolution in Sectional Matrix Systems: The Quad Matrix System™

Spring Lake, MI – Garrison Dental Solutions, a leader in advanced matrix systems, is excited to announce the launch of the Quad Matrix System, an innovative solution designed to empower clinicians in tackling the most challenging Class II restorations.

Beyond Conventional Class II Restorations

The Quad Matrix System overcomes the limitations of traditional matrix systems, effectively addressing complex scenarios such as back-to-back restorations, deep caries, and multi-tooth restorations with precision and efficiency. Whether it's a single restoration or an entire quadrant of Class IIs, the Quad system provides the versatility and performance clinicians need.

Clinician-Inspired, Garrison-Engineered

Conceived by Dr. Matt Burton, a practicing dentist, and meticulously brought to life by Garrison's engineering expertise, the Quad Matrix System is backed by thousands of hours of clinical validation. It addresses common challenges with unique features that enhance tooth separation, matrix stabilization, and cervical sealing.

Redefining Performance with an Advanced Design

- **Driver-Tip Quad Rings:** The asymmetrical tip design maximizes ring adaptation to effectively address complex restorations.
- **Split-Tip Quad Wedges:** Designed to work synergistically with the Quad ring to reinforce gingival seal, and support a natural emergence profile, creating superior contact areas.
- **Firm and Dead-Soft Bands:** Comprehensive range of firm and dead-soft bands to accommodate a variety of clinical scenarios, from single Class IIs to back-to-backs, ensuring superior results across all procedures.

Comprehensive Solutions for Complex Restorations

The Quad Matrix System supports a wide range of clinical scenarios:

- **Back-to-Back Restorations:** Simultaneously restore both sides of adjacent Class IIs using two matrix bands, one Quad Wedge, and one Quad Ring.
- **Deeper Caries:** Achieve a secure seal for subgingival restorations with Firm Bands and Quad Wedges enhanced by the driver-tip design.
- **Multi-Tooth Restorations:** Use Universal (blue) and Tall (orange) Quad Rings to create ideal contours and contacts.

Setting a New Standard in Dental Restorations

The Quad Matrix System by Garrison Dental significantly advances how clinicians approach Class II restorations. Its innovative features, supported by a fully integrated design, provide unmatched versatility, efficiency, and precision for even the most complex cases.

Availability

The Quad Matrix System is now available. For more information, visit www.garrisondental.com or contact your Garrison representative.

About Garrison Dental Solutions

Garrison Dental Solutions has been a pioneer in dental matrix systems for over 25 years. Dedicated to innovation and clinician success, Garrison develops products that enhance patient care and streamline workflow.